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- Lieu, J. E. (2004). Speech-language and educational consequences of unilateral hearing loss in children. *Archives of Otolaryngology - Head and Neck Surgery*, 130(5), 524–530.
- Mülder, H. E., & Smaka, C. (2013). Interview with Dr. Hans E. Mülder, Director Marketing and Senior Audiologist at Phonak Communications, Phonak Headquarters, Switzerland. *AudiologyOnline* (Online), retrieved from <http://www.audiologyonline.com/interviews/interviewwithdrs-hans-e-11727>, accessed June 17th, 2019. Phonak AG. (2013).
- Phonak Insight | Roger Pen—Bridging the understanding gap. Retrieved from https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/evidence/white_paper/documents/technical_paper/Insight_Roger_Pen_028-0933.pdf, accessed June 17th, 2019.
- Phonak AG. (2014). Fitting Guide Roger and Cochlear sound processors Nucleus 5 and Nucleus 6. (Online). Retrieved from https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/products_solutions/wireless_accessories/roger_receivers/documents/fitting_guide_roger_cochlear_nucleus%205_6.pdf, accessed June 17th, 2019.
- Schafer, E. C., Thibodeau, L. M. (2006). Speech recognition in noise in children with cochlear implants while listening in bilateral, bimodal and FM system arrangements. *American Journal of Audiology*, 15(2), 114–126.
- Távora-Vieira, D., De Ceulaer, G., Govaerts, P. J., Rajan, G. P. (2015). Cochlear implantation improves localization ability in patients with unilateral deafness. *Ear & Hearing*, 36(3), 93–98.
- Thibodeau, L. (2014). Comparison of speech recognition with adaptive digital and FM remote microphone hearing assistance technology by listeners who use hearing aids. *American Journal of Audiology*, 23(2), 201–210.
- Vermeire, K., Van de Heyning, P. (2009). Binaural hearing after cochlear implantation in subjects with unilateral sensorineural deafness and tinnitus. *Audiology & Neuro-Otology*, 14(3), 163–171.
- Wagener, K., Kühnel, V., & Kollmeier, B. (1999a). Entwicklung und Evaluation eines Satztests in deutscher Sprache I: design des Oldenburger Satztests. *Z Audiology*, 38(1), 4–15.
- Wagener, K., Brand, T., & Kollmeier, B. (1999b). Entwicklung und Evaluation eines Satztests in deutscher Sprache III: evaluation des Oldenburger Satztests. *Z Audiology*, 38(3), 86–95.
- Wesarg, T. Arndt, S., Wiebe, K., Schmid, F., Huber, A., Hans, E., Mülder, H., Laszig, R., Aschendorff, A., & Speck, I. (2018). Speech recognition in noise in single-sided deaf cochlear implant recipients using digital remote wireless microphone technology. *Journal of the American Academy of Audiology*, DOI: jaaa17131q.
- Wie, O. B., Pripp, A. H., & Tvette, O. (2010). Unilateral deafness in adults: effects on communication and social interaction. *Annals of Otology, Rhinology & Laryngology*, 119(11), 772–781.
- Wolfe, J., Morais, M., Schafer, E., Mills, E., Mülder, H. E., Goldbeck, F., Marquis, F., John, A., Hudson, M., Peters, B. R., & Lianos, L. (2013). Evaluation of speech recognition of cochlear implant recipients using a personal digital adaptive radio frequency system. *Journal of the American Academy of Audiology*, 24(8), 714–724.

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