Pediatric Focus 1

Январь 2017 – Dawna Lewis и Marlene Bagatto

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McCreery, R., Venediktov, R., Coleman, J., & Leech, H. (2012). An evidence-based systematic review of directional microphones and digital noise reduction hearing aids in school-age children with hearing loss. *American Journal of Audiology*, 21, 295–312

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American Academy of Audiology. (2013). Clinical practice guidelines on pediatric amplification. Retrieved from www.audiology.org/resources/documentlibrary/Documents/PediatricAmplificationGuidelines.pdf

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Protocol for the Provision of Amplification within the Ontario Infant Hearing Program (MCYS, 2014; <http://www.dsl.io/wp-content/uploads/2015/05/DSL5-Pediatric-Protocol.2014.01.pdf>)

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Anderson, K. (2001). Kids in noisy classrooms: What does the research really say? *Journal of Educational Audiology*, 9, 21-33.

Bradley, J., & Sato, H. (2008). The intelligibility of speech in elementary school classrooms. *Journal of the Acoustical Society of America*, 123, 2078–2086.

Dockrell, J. & Shield, B. (2006). Acoustical barriers in classrooms; the impact of noise on performance in the classroom. *British Educational Research Journal*, 32, 509-525.

Elliott, L. (1979). Performance of children aged 9 to 17 years on a test of speech intelligibility in noise using sentence material with controlled word predictability. *Journal of the Acoustical Society of America*, 66, 651–653.

Fallon, M., Trehub, S., & Schneider, B. (2000). Children's perception of speech in multitalker babble. *Journal of the Acoustical Society of America*, 108, 3023–3029.

Fallon, M., Trehub, S., & Schneider, B. (2002). Children's use of semantic cues in degraded listening environments. *Journal of the Acoustical Society of America*, 111(5 Pt 1), 2242–2249.

Hygge, S., Evans, G., & Monika, B. (2002). A prospective study of some effects of aircraft noise on cognitive performance in school children. *Psychological Science*, 13, 469-474.

Jamieson, D. G., Kranjc, G., Yu, K., et al. (2004). Speech intelligibility of young school-aged children in the presence of real-life classroom noise. *Journal of the American Academy of Audiology*, 15, 508–517.

Johnson, C. (2000). Children's phoneme identification in reverberation and noise. *Journal of Speech, Language, and Hearing Research*, 43, 144–157.

Klatte, M., Hellbruck, J., Seidel, J., & Leistner, P. (2010). Effects of classroom acoustics on performance and well-being in elementary school children. *Environment & Behavior*, 42, 659-692.

Klatte, M., Lachmann, T., & Meis, M. (2010). Effects of noise and reverberation on speech perception and listening comprehension of children and adults in a classroom-like setting. *Noise & Health*, 12, 270-282.

McCreery, R., Ito, R., Spratford, M., Lewis, D., Hoover, B., & Stelmachowicz, P. (2010). Performance–intensity functions for normal-hearing adults and children using computer-aided speech perception assessment. *Ear & Hearing*, 31, 95-101.

McCreery, R. & Stelmachowicz, P. (2013). The effects of limited bandwidth and noise on verbal processing time and word recall in normal-hearing children. *Ear & Hearing*, 34, 585-591.

Neuman, A., & Hochberg, I. (1983). Children's perception of speech in reverberation. *Journal of the Acoustical Society of America*, 73, 2145-2149.

Neuman, A., Wroblewski, M., Hajicek, J., et al. (2010). Combined effects of noise and reverberation on speech recognition performance of normal hearing children and adults. *Ear & Hearing*, 31, 336-344.

Nittrouer, S. & Boothroyd, S. (1990). Context effects in phoneme and word recognition by young children and older adults. *Journal of the Acoustical Society of America*, 87, 2705-2715.

Shield, B., & Dockrell, J. (2008). The effects of environmental and classroom noise on the academic attainments of primary school children. *Journal of the Acoustical Society of America*, 123, 44-133. doi:10.1121/1.2812596.

Wroblewski, M., Lewis, D., Valente, D., et al. (2012). Effects of reverberation on speech recognition in stationary and modulated noise by school aged children and young adults. *Ear & Hearing*, 33, 731-744.

Yang, W., & Bradley, J. S. (2009). Effects of room acoustics on the intelligibility of speech in classrooms for young children. *Journal of the Acoustical Society of America*, 125, 922-933.

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Ching, T. & Hill, M. (2007). The Parents' Evaluation of Aural/Oral Performance of Children (PEACH) scale: normative data. *Journal of the American Academy of Audiology*, 18, 220-235.

Crandell, C. & Smaldino, J. (1994). An update of classroom acoustics for children with hearing impairment. *Volta Review*, 96, 291-306.

Finitzo-Hieber, T. & Tillman, T. (1978). Room acoustics effects on monosyllabic word discrimination ability for normal and hearing-impaired children. *Journal of the Acoustical Society of America*, 21, 440-458.

Hicks, C.B. & Tharpe, A.M. (2002). Listening effort and fatigue in school-age children with and without hearing loss. *Journal of Speech, Language, and Hearing Research*, 45, 573-584.

Iglehart, F. (2009). Combined effects of classroom reverberation and noise on speech perception by students with typical and impaired hearing. Paper presented at Inter-Noise 2009, Ottawa, Ontario, Canada.

Leibold L., Hillock-Dunn A., Duncan N., Roush P., Buss E. (2013). Influence of hearing loss on children's identification of spondee words in a speech-shaped noise or a two-talker masker. *Ear & Hearing*, 34, 575-584.

Stelmachowicz, P., Pittman, A., Hoover, B., & Lewis, D. (2001). Effect of stimulus bandwidth on the perception of /s/ in normal- and hearing-impaired children and adults. *Journal of the Acoustical Society of America*, 110, 2183-2190.

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Anderson, K., & Goldstein, H. (2004). Speech perception benefits of FM and infrared devices to children with hearing aids in a typical classroom. *Language, Speech, and Hearing Services in Schools*, 35, 169-184.

Anderson, K., Goldstein, H., Colodzin, L., & Iglehart, F. (2005). Benefit of S/N enhancing devices to speech perception of children listening in a typical classroom with hearing aids or a cochlear implant. *Journal of Educational Audiology*, 12, 14-28.

Hawkins, D. (1984). Comparisons of speech recognition in noise by mildly-to-moderately hearing-impaired children using hearing aids and FM systems. *Journal of Speech and Hearing Disorders*, 49, 409-418.

Pittman, A., Lewis D., Hoover B., Stelmachowicz, P. (1999). Recognition performance for four combinations of FM system and hearing aid microphone signals in adverse listening conditions. *Ear & Hearing*, 20, 279-289.

Thibodeau, L. (2010). Benefits of adaptive FM systems on speech recognition in noise for listeners who use hearing aids. *American Journal of Audiology*, 19, 36-45.

Wolfe, J., Morais, M., Neumann, S., et al. (2013). Evaluation of Speech Recognition with Personal FM and Classroom Audio Distribution Systems. *Journal of Educational Audiology*, 19, 65-79.

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Bentler, R. (2005). Effectiveness of directional microphones and noise reduction schemes in hearing aids: a systematic review of the evidence. *Journal of the American Academy of Audiology*, 16, 473-484.

Kuk, F., Keenan, D., Lau C., & Ludvigsen, C. (2005). Performance of a fully adaptive directional microphone to signals presented from various azimuths. *Journal of the American Academy of Audiology*, 16, 333-347.

Mackenzie, E. & Lutman, M. (2005). Speech recognition and comfort using hearing instruments with adaptive directional characteristics in asymmetric listening conditions. *Ear & Hearing*, 26, 669-679.

Mueller, H., Weber, J., & Hornsby, B. (2006). The effects of digital noise reduction on the acceptance of background noise. *Trends in Amplification*, 10, 83-93.

Palmer, C., Bentler, R., & Mueller, H. (2006). Amplification with digital noise reduction and the perception of annoying and aversive sounds. *Trends in Amplification*, 10, 95-104.

Ricketts, T., Henry, P., & Gnewikow, D. (2003) Full time directional versus user selectable microphone modes in hearing aids. *Ear & Hearing*, 24, 424-39.

Ricketts, T., Henry, P., & Hornsby, B. (2005). Application of frequency importance functions to directivity for prediction of benefit in uniform fields. *Ear & Hearing*, 26, 473-486.

Ricketts, T. & Hornsby, B. (2005). Sound quality measures for speech in noise through a commercial hearing aid implementing digital noise reduction. *Journal of the American Academy of Audiology*, 16, 270-277.

Ricketts, T. & Hornsby, B. (2006). Directional hearing aid benefit in listeners with severe hearing loss. *International Journal of Audiology*, 45, 190-197.

Walden, B., Surr, R., Cord, M., et al. (2007). The robustness of hearing aid microphone preferences in everyday listening environments. *Journal of the American Academy of Audiology*, 18, 358-379.

Yuen, K., Kam, A., & Lau, P. (2006). Comparative performance of an adaptive directional microphone system and a multichannel noise reduction system. *Journal of the American Academy of Audiology*, 17, 241-252.

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American Academy of Audiology. (2013). Clinical practice guidelines on pediatric amplification. Retrieved from www.audiology.org/resources/documentlibrary/Documents/PediatricAmplificationGuidelines.pdf.

Aurimmo, J., Kuk, F., Lau, C., et al. (2009). Efficacy of an adaptive directional microphone and a noise reduction system for school aged children. *Journal of Educational Audiology*, 15, 15–27.

Gravel, J. S., Fausel, N., Liskow, C., & Chobot, J. (1999). Children's speech recognition in noise using omnidirectional and dual-microphone hearing aid technology. *Ear & Hearing*, 20, 1–11.

Hawkins, D. B. (1984). Comparisons of speech recognition in noise by mildly-to-moderately hearing-impaired children using hearing aids and FM systems. *Journal of Speech and Hearing Disorders*, 49, 409–418.

Kuk, F., Kollofski, C., Brown, S., Melum, A., & Rosenthal, A. (1999). Use of a digital hearing aid with directional microphones in school-aged children. *Journal of the American Academy of Audiology*, 10, 535–548.

McCreery, R., Venediktov, R., Coleman, J., & Leech, H. (2012). An evidence-based systematic review of directional microphones and digital noise reduction hearing aids in school-age children with hearing loss. *American Journal of Audiology*, 21, 295–312.

Pittman, A. & Hiipakka, M. (2013). Hearing impaired children's preference for, and performance with, four combinations of directional microphone and digital noise reduction technology. *Journal of the American Academy of Audiology*, 24, 832-44.

Ricketts, T.A. & Galster, J. (2008). Head angle and elevation in classroom environments: implications for amplification. *Journal of Speech, Language, and Hearing Research*, 51, 516-525.

Ricketts, T., Galster, J., & Tharpe, A. M. (2007). Directional benefit in simulated classroom environments. *American Journal of Audiology*, 16, 130–143.

Ricketts, T. & Picou, E. (2013). Speech recognition for bilaterally asymmetric and symmetric hearing aid microphone modes in simulated classroom environments. *Ear & Hearing*, 34, 601-609.

Wouters, J., Litière, L., & van Wieringen, A. (1999). Speech intelligibility in noisy environments with one- and two-microphone hearing aids. *Audiology*, 38, 91–98.

Protocol for the Provision of Amplification (infants and preschool children; Ontario IHP, 2014). Retrieved from <http://www.dslio.com/wp-content/uploads/2015/05/DSL5-Pediatric-Protocol.2014.01.pdf>

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