

Rediscover the wonder of sound



Phonak CROS P

A Paradise solution for unilateral hearing loss

A Sonova brand

PHONAK
life is on

Rediscover
the wonder
of sound
around you
with Phonak
CROS P

Well-Hearing is Well-Being

Hearing with one ear is challenging. If you're in a noisy room and the person you're conversing with is positioned on the side of your poorer hearing ear, it can be difficult to follow the conversation or even acknowledge that someone is talking. This is where Phonak CROS P, a solution for unilateral hearing loss, can really make a difference to the quality of your hearing and your overall well-being. Because at Phonak we truly believe that hearing well fosters easier engagement and stronger connections.

Benefits of CROS devices for unilateral hearing loss:^{1*}

- Follow** conversations wherever they come from.
- Hear speech better from your poorer side even when in noise
- Gain greater environmental awareness of sounds coming into your poorer hearing ear



¹ Stewart, E. & Woodward, J. (2021). Out of the (Head) Shadow: A Systematic Review of CROS/BiCROS Literature. *The Hearing Review*, 28(8),22-25. Retrieved from <https://www.hearingreview.com/inside-hearing/research/head>, accessed August 9, 2021.

* When using a CROS system compared to unaided.

** Follow - means awareness of conversation for the wearer to attend to.

Make the good even better!

CROS P will transfer sound and voices, enabling your better hearing ear to access sounds from both sides.

If you have normal hearing or aidable hearing loss in one ear and no usable hearing in the other ear, CROS P is designed for you. You'll be able to hear your partner or colleague in a noisy cafe even though they are talking on your poorer side.

CROS stands for 'Contralateral Routing of Signal'

It consists of two parts:

- ① The CROS with a microphone that picks up sounds and voices presented to the poorer hearing side and wirelessly transmits them to the hearing aid.
- ② The Audéo Paradise hearing aid that receives the sounds from the poorer hearing side and plays it to your better hearing ear.

Wireless transfer



① Phonak CROS P



② Phonak Audéo P hearing aid





A Paradise solution for unilateral hearing loss

Phonak CROS P is based on the proven performance of Phonak Paradise technology. Paradise delivers an unrivaled hearing experience* in even more situations. When CROS P is paired with an Audéo Paradise hearing aid, it enables you to hear sounds around you and follow*** conversations, regardless of the direction they come from.^{1**}

CROS P utilizes AutoSense OS™ 4.0 – the proprietary operating system from Phonak. AutoSense OS automatically senses your surroundings and adjusts your hearing aid accordingly. It orchestrates various features such as Binaural VoiceStream Technology™, enabling CROS P to wirelessly transfer sounds from the poorer hearing side to the hearing aid on the better hearing ear. This combination of Phonak features allows you to hear your loved ones regardless if they are talking to your poorer hearing side so you can be in the center of the conversation again.

* Compared to Phonak Audéo™ Marvel.

** When using a CROS system compared to unaided.

*** Follow - means awareness of conversation for the wearer to attend to.

¹ Stewart, E. & Woodward, J. (2021). Out of the (Head) Shadow: A Systematic Review of CROS/BiCROS Literature. *The Hearing Review*, 28(8), 22-25. Retrieved from <https://www.hearingreview.com/inside-hearing/research/head>, accessed August 9, 2021.



Connectivity is key

Your CROS P system has universal connectivity to smartphones, TV, Roger™ microphones and more. Allowing you to enjoy hands-free calls or stream media content from iOS®, Android™ and other Bluetooth® enabled devices. This solution is easy to use and includes a rechargeable model. Just turn it on and the rest is automatic.

Highlights of Phonak CROS P:

- Proven Paradise performance
- Improved speech understanding in noisy environments^{1*}
- Follow** conversations wherever they are coming from^{1*}
- Universal connectivity to Bluetooth enabled smartphones, TV, Roger microphones and more
- Easy to use
- Rechargeable



iOS® is a registered trademark or trademark of Cisco Systems, Inc. and/or its affiliates in the United States and certain other countries.

Android is a trademark of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Sonova AG is under license.

¹ Stewart, E. & Woodward, J. (2021). Out of the (Head) Shadow: A Systematic Review of CROS/BiCROS Literature. *The Hearing Review*, 28(8), 22-25. Retrieved from <https://www.hearingreview.com/inside-hearing/research/head>, accessed August 9, 2021.

* When using a CROS system compared to unaided.

** Follow - means awareness of conversation for the wearer to attend to.

Phonak CROS P in different situations

Conversations in noisy places

CROS P helps you to hear better in noisy environments¹ especially when speech is coming to your poorer hearing side. Dedicated features in Paradise help you to **focus on speech when you need it most.**

Taking phone calls with a tap

Your CROS P system can connect to more than 6800 smartphone models² plus Bluetooth® enabled devices, turning your hearing solution into wireless headphones. **By simply tapping your ear, you can even take phone calls, pause and resume streamed music and activate voice assistant apps.**

Walking down the street

CROS P will keep you aware of sounds when you're walking down the street because it **enables you to hear sounds coming from both directions.** This means you can **be in the center of the conversation** without having to consider the best position to stand in. CROS P enables you to hear voices wherever they are coming from.

¹ Stewart, E. & Woodward, J. (2021). Out of the (Head) Shadow: A Systematic Review of CROS/BiCROS Literature. *The Hearing Review*, 28(8), 22-25. Retrieved from <https://www.hearingreview.com/inside-hearing/research/head>, accessed August 9, 2021.

² [phonak.com/compatibility](https://www.phonak.com/compatibility), over 6800 compatible models, March 2021



The right or left fit!

CROS P is available in two discreet models that fit comfortably behind your ear.

- CROS P-13 offers extended battery life with a replaceable zinc-air battery.
- CROS P-R offers 12.5 hours of hearing from both sides and Tap Control with no battery hassles.



Phonak CROS P-13



Phonak CROS P-R

For your good ear, you can pair with an Audéo P-13, P-R or P-RT hearing aid.

Color coordinated

CROS P is available in eight colors to suit your needs:



Sand Beige
P1



Sandalwood
P3



Chestnut
P4



Champagne
P5



Silver Gray
P6



Graphite Gray
P7



Velvet Black
P8



Beige
01

Get the most out of your Phonak CROS P

Whether you're getting a CROS for the very first time or upgrading from a previous CROS system, here are our top tips.

① Talk to your hearing care professional (HCP)

Your HCP will ensure that your settings are comfortable. If there are any situations where you feel the sound was not quite right – make a note and go back to your HCP to get CROS P set up in a way that best fits you.

② Give it time

Adjusting to anything new takes time; whether it's a new pair of shoes or a new car. Adjusting to CROS P is no different. Everything will sound different at the beginning so be patient with yourself.

③ Try different environments

Try out familiar situations and then start trying out new environments. Is there somewhere you've avoided because of your hearing loss? See how CROS P makes a difference to your comfort level and enjoyment.

④ **Tell your friends and family**

Let people know that you're adjusting to new hearing aids and it may take some time to get used to.

⑤ **Give yourself a break**

If you've had unilateral hearing loss for a long time, it may take your brain some time to adjust. That is normal. Give yourself a break when needed but remember that the more you use your new CROS P, the better it will be.

life is on

At Phonak, we believe that well-hearing equates to well-being and is essential to living life to the fullest. For more than 70 years, we have remained passionate about creating a world where 'life is on' for everyone. Our innovative hearing solutions are designed for people of all ages and all degrees of hearing loss, to connect socially, thrive mentally and emotionally.

www.phonak.com



7 613389 520262

028-2296-02/V2.00/2021-09/Jt © 2021 Sonova AG. All rights reserved.



sonova
HEAR THE WORLD



Sonova AG · Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland



Sonova Deutschland GmbH · Max-Eyth-Str. 20
70736 Fellbach-Oeffingen · Germany