

Phonak

Life Charger™

Technical Data

Dimensions (L x W x H)	94 x 51 x 62 mm (3.7 x 2 x 2.4")
Weight	145 g (5.1 oz)
Color	Champagne
Operating conditions	5° to +30° Celsius (41° to 86° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation	Up to 70% (non condensing)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 85% (non condensing)
Atmospheric pressure during use	800 hPa to 1060 hPa
Atmospheric pressure during storage and transportation	500 hPa to 1060 hPa
Included Parts	Power supply USB-C cable (for charging only)
External Interfaces	USB-C port
Input	Type: DC, stabilized 5 V +/- 10% Max. current: 1 A



Sonova AG · Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland
www.phonak.com

A Sonova brand

PHONAK
life is on