

Initial diagnostic testing of infants who fail newborn screening

Martyn Hyde

Professor, Otolaryngology/S-LP, University of Toronto
Assoc Director, Hearing, Balance & Speech Dept, Mt Sinai Hosp
Consultant, Ontario Infant Hearing Program &
British Columbia Early Hearing Program

mhyde@mtsinai.on.ca

1

Best textbook chapter ever written on this subject !

Thanks to:

Dr David Stapells
University of British Columbia
Vancouver, Canada

Frequency-Specific Threshold Assessment in Young Infants
Using the Transient ABR and the Brainstem ASSR

Chapter 21, in: RC Seewald & AM Tharpe (Eds)
Comprehensive Handbook of Pediatric Audiology
San Diego, Plural Publishing

2

Challenge: after AABR screen fail or direct high-risk referral:

- Significant hearing loss? YES / NO
- If YES, what type?
 - Conductive (T or P) Conventional cochlear ANSD Brainstem neural
 - Any mixture of these !
- If NOT Conductive (T), ANSD or neural:
Thresholds sufficient to specify amplification, Implant candidacy, baseline for progression

3

Additional challenges

- Validity & accuracy, given limited data in young infants with proven permanent hearing loss, especially for ASSR & for bone-conduction
- Efficiency & completeness, given three pressures: access, baby EEG state & early intervention

4

Present diagnostic tools

- Tonepip ABR (air & bone conduction)
- Brainstem (80 Hz) ASSR: single/multiple-frequency
- Click ABR & cochlear microphonic (CM)
- (Long-latency cortical potentials)

- OAE (Distortion Product or Transient): 1-4 kHz
- Tympanometry: 1 kHz probe < 6 months
- Middle-ear muscle reflexes: ipsi, 1k, wide-band



5


AABR screen will NOT fail:

- ANY abnormal audiogram with ANY threshold better than 30 dBHL in the range 1-8 kHz


- These include:
 - Sloping losses normal at 1kHz or above
 - Steep high-frequency losses above 1kHz
 - Reverse-slope losses normal at 8kHz or below
 - U-shaped losses normal at 1kHz or 8kHz

6






Click ABR present at 30 dBnHL




Click ABR threshold 60 dBnHL
Same amplification...???



11

Click ABR threshold 50 dBnHL
Monitoring: no progression..???




12

Clicks are an otoneurologic tool, not a sufficient threshold tool

- Would you do an adult hearing test with white noise only?
- Infants deserve the best possible audiometry
- Clicks have many technical deficiencies & may soon be replaced by chirps

13

Best published data to date + long experience




Stapells, Gravel & Martin, 1995

Estimated Hearing Level (eHL)

- Determine tone-ABR threshold:
ABR present at x dBnHL, absent at x-10
- Correction factor C is normative median difference
(ABR threshold – HL threshold)
- eHL threshold = ABR threshold – C
- AC 500 15 dB 1k 10 dB 2k 5 dB 4k 0 dB

15

Minimum intensity levels (25 dB HL) Recording time < 10 mins



16

Successful diagnostic ABR testing

- Sufficient caseload to develop/maintain skills
- 4-8 weeks corrected age, 4 weeks after screen
- Natural sleep, arrive tired & hungry, feed after electrode attachment
- TWO insert transducers, whenever possible
- Hand-held BC if trained, otherwise use Velcro band
- High-efficiency (top-down, progressive) strategy:
Every test condition is chosen so its result will make the greatest difference to baby MANAGEMENT if the test stopped at that point...

17

Minimum test for AABR failures?


Loss probability: No risk 0.1, High risk 0.4

- Statistical false positive AABR – N hearing
- ABR detection failure – high EEG noise levels
- Transient hearing loss/debris/probe blockage
- Permanent hearing loss (1-8 kHz)/ANSD
- Test options: 2k, 4k, 1k, 500
2k, 500 or 4k, 1k
2k only !?
- ABR at minimum required level? Yes: STOP
No: ENTER EFFICIENT AUDIOGRAM STRATEGY

18

Minimum test for direct Dx referral (eg very high risk, no AABR screen)

- ANY hearing loss profile is possible, so:
- ENTER EFFICIENT AUDIOGRAM STRATEGY



19


Efficient loss type strategy

- Try OAE first if baby not asleep (esp. high risk)
- If OAE 2k present & 2k ABR absent at 90 dBnHL, ANSD is present
- If first ear 2k 'normal', change ears immediately
- If first ear 2k threshold ≥ 40 dBnHL go to BC
- If BC 2k ABR present at min level, stop


Efficient threshold strategy

- Start with 2kHz at minimum level & go up in LARGE steps (20-30 dB)
Never do input-output function with 10 dB steps !
- Use 1-2000 sweeps with repetition if needed
- Bracket threshold with 10 dB step, repeated averages & use residual noise level (RNL)
- Aim to locate threshold in ≤ 6 averages


Example of inefficiency




False +ve & low-noise -ve



BC: Inference of responding cochlea




2 kHz BC ipsi-contra in a unilateral



Common error in ABR detection

- There are **THREE** decision outcomes, not **TWO**: Present, Absent, and **INDETERMINATE**
- If EEG noise level is not low enough, you cannot and must not say that response is absent or present (different criteria for each)
- Guessing is much worse than saying 'I don't know'
- Should use residual noise level if available

Not so easy




Using SNR & residual noise (IHS)

Present




SNR 1.26 RN 0.07 μ V

Absent



SNR 0.64 RN 0.04 μ V

Indeterminate



SNR 0.16 RN 0.15 μ V

ASSR: some comments


- 80Hz ASSR is probably ABR V-V' – difficult to prove
- High stimulus rate is efficient, objective detection is convenient, immunity to 60 Hz interference is useful
- 'Better frequency-specificity' NOT proven in infants with definite SNHL
- 'Can measure greater hearing losses' is probably an artifact of dBnHL vs HL (or vestibular response)

Tone-ABR vs s-ASSR vs m-ASSR

- m-ASSR less efficient for frequency-intensity strategy
- In young infants with proven SNHL:
 - Very few *good* studies of s- or m-ASSR accuracy
 - No normative studies of BC
 - No good studies of relative efficiency
 - Insufficient data on stimulus interactions
 - Optimal test parameters are not well-established

Behavioural (VRA) vs m-ASSR, 21 infants

Luts et al *Audiol Neurotol* 2006;11:24-37




BUT...we *can* say that in infants:

- Multiple ASSR (one or two ears, AC) is a valid & efficient tool to prove absence of significant hearing loss, given response at:

	500 Hz	1 kHz	2kHz	4kHz
dBHL	50	45	40	40

- Recent data show a strong relationship between multiple-ASSR & tone-ABR thresholds (Van Maanen & Stapells, *JAAA* 2010)

Multiple ASSR (AC) vs tone-ABR in infants with hearing loss (*Van Maanen & Stapells 2010*)



Mensajes para llevar...

- Clicks are not good enough
- Tone ABR AC & BC is valid & practicable
High-efficiency standard protocol is essential
- DO NOT need local norms. Google: 'BC Early Hearing Program Diagnostic Audiology Protocol'
- ASSR is efficient for normal/abnormal decision
- Insufficient ASSR data for complete audiometry

