

# Phonak CROS

## Product information



Phonak CROS is a CROS/BiCROS solution for people with unaidable hearing loss on one side. It consists only of two parts – the Phonak CROS transmitter and a wireless Phonak Quest or Spice+ hearing instrument. The Phonak CROS transmitter includes two BTE models and different ITE form factors, each with the choice of a size 13 or 312 battery. They can be combined with all form factors in the premium, advanced and standard performance level, depending on individual needs and preferences.

For more information visit [www.phonakpro.com](http://www.phonakpro.com)

# Product description



Phonak CROS



Phonak CROS H20



Phonak CROS 312



Phonak CROS 13

Frequency range (Hz)	2cc	Audible frequency range depending on hearing instrument	Audible frequency range depending on hearing instrument	Audible frequency range depending on hearing instrument	Audible frequency range depending on hearing instrument
Battery size		312	13	312	13
Working current (mA)		~3,8	~3,8	~3,5	~3,5
Shell style				ITC Half shell Full shell	ITC Half shell Full shell
Microphone technology		Real Ear Sound or omnidirectional	Real Ear Sound or omnidirectional	omnidirectional	omnidirectional
Push button (Program / Volume)		•	•	•	•
Volume control					
Telecoil					
Nano coating		•	•		
WaterResistant – IP67 <sup>1</sup>			•		
Direct Audio Input (DAI)					
Secured battery compartment			•		

<sup>1</sup> IP67 indicates that the hearing aid was not damaged beyond repair after 8 hours in a dust chamber, nor after being immersed in 1 meter (3 feet) of water for 30 minutes as defined by IEC60529

# Housing colors

## Phonak Quest color range

Phonak CROS / Phonak CROS H2O

### Hair & Skin

<b>P1</b> Sand Beige	<b>P2</b> Amber Beige	<b>P3</b> Sandalwood	<b>P4</b> Chestnut	<b>P5</b> Champagne	<b>P6</b> Silver Gray	<b>P7</b> Graphite Gray	<b>P8</b> Velvet Black
-------------------------	--------------------------	-------------------------	-----------------------	------------------------	--------------------------	----------------------------	---------------------------

### Fashion

<b>P9</b> Ruby
-------------------

<b>Q1</b> Petrol
---------------------

### Fun

<b>13</b> Pure Transparent
----------------------------------

<b>38</b> Purple Transparent
------------------------------------

<b>37</b> Blue Transparent
----------------------------------

### Traditional

<b>01</b> Beige
--------------------

## Phonak Spice+ color range

### Phonak CROS

<b>F6</b> Beach Beige	<b>F7</b> Mocha Taupe	<b>G1</b> Mid Cloud Gray	<b>F8</b> Chocolate Brown	<b>F9</b> Midnight Black	<b>C0</b> Silver Shadow	<b>G3</b> Ocean Blue
--------------------------	--------------------------	--------------------------------	---------------------------------	-----------------------------	----------------------------	-------------------------

### Phonak CROS H2O

<b>01</b> Beige	<b>F7</b> Mocha Taupe	<b>G1</b> Mid Cloud Gray	<b>E5</b> Chocolate	<b>06</b> Black	<b>70</b> Palladium Black	<b>G3</b> Ocean Blue
--------------------	--------------------------	--------------------------------	------------------------	--------------------	---------------------------------	-------------------------

## Faceplate colors

### Phonak CROS 312 / Phonak CROS 13

<b>14</b> Tan	<b>22</b> Cocoa	<b>26</b> Pink	<b>28</b> Brown
------------------	--------------------	-------------------	--------------------

## Shell colors





### Phonak CROS 312 / Phonak CROS 13

<b>14</b> Tan	<b>22</b> Cocoa	<b>26</b> Pink	<b>28</b> Brown	<b>19</b> White	<b>21</b> Clear Transparent	<b>36</b> Red Transparent	<b>37</b> Blue Transparent
------------------	--------------------	-------------------	--------------------	--------------------	-----------------------------------	---------------------------------	----------------------------------

# Acoustic coupling

Universal retention	Phonak CROS Retention	Available as left and right
	Phonak CROS SlimTube	Available in left, right and four different lengths ( 0, 1, 2, 3)
Customized retention	Phonak CROS Tip	Custom made earpiece available in hard and soft material
Earpiece	Dome	Open dome available in three sizes (S, M, L)
		Closed dome available in three sizes (S, M, L)
		Power dome available in three sizes (S, M, L)

# Wireless communication compatibility

Digital Wireless Accessories <sup>1</sup>	Phonak CROS	Phonak CROS H20	Phonak CROS 312	Phonak CROS 13
 Phonak PilotOne	•	•	•	•
 Phonak ComPilot	•	•	•	•
 Phonak TVLink S	•	•	•	•
 Phonak RemoteMic	• <sup>2</sup>	• <sup>2</sup>	• <sup>2</sup>	• <sup>2</sup>

<sup>1</sup> **Attention:** When wireless accessories or FM systems are used (excluding remote controls), sound from the CROS transmitter side will not be heard in the hearing instrument.

<sup>2</sup> Requires ComPilot

# Fitting

Software	Phonak Target™ 2.1 or higher
Interfaces	iCube, NOAHLink, HI-PRO, HI-PRO 2

# Key features

	Phonak CROS	Phonak CROS H20	Phonak CROS 312	Phonak CROS 13
Binaural VoiceStream Technology™	•	•	•	•
Real Ear Sound	•	•		
QuickSync	•	•	•	•
Warning signal				
CROS connected	•	•	•	•
CROS disconnected	•	•	•	•
Battery warning	•	•	•	•
Wireless fitting	•	•	•	•

# Hearing instrument compatibility matrix

	Phonak Quest (wireless hearing instruments)				Phonak Spice+ (wireless hearing instruments)		
	Phonak Bolero Q	Phonak Virto Q	Phonak Naída Q <sup>1</sup>	Phonak Audéo Q <sup>1</sup>	BTE <sup>2</sup>	ITE <sup>2</sup>	Audéo S <sup>3</sup>
Phonak CROS	M312				microM		
Phonak CROS H20	M13	312	RIC	312	M H20	312	
Phonak CROS 312	P	13	SP	312T	microP	312 UZ	SMART
Phonak CROS 13	SP		UP		SP	13 UZ	YES
					UP		

<sup>1</sup> Not compatible with Q30

<sup>2</sup> Not compatible with Phonak Dalia

<sup>3</sup> Not compatible with Audéo SMART I

---

# Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence.  
Live without limit. Life is on.

[www.phonakpro.com](http://www.phonakpro.com)