

Quick reference guide

		Technical Information						Phonak AccessLine						
		Battery	Ear Simulator		2 cc		Acoustic Coupling	Categories of Hearing Impairment	Accessories			FM Receiver		FM Transmitter
			Max. Gain	Max. MPO	Max. Gain	Max. MPO	SlimTube HE Hook SlimTip/SlimTip soft Open/closed/power dome	Mild hearing loss Moderate hearing loss Severe hearing loss Profound hearing loss	Phonak PilotOne myPilot iCom Click'nTalk Phonak TVLink iCube	MyLink+	MLiS or MLi MLi3i AS13	inspiro SmartLink+ ZoomLink+ EasyLink+		
Phonak Ambra Phonak Solana Phonak Cassia	Petite	10	51	120	45	118								
	microM	312	59	128	54	126								
	microP	13	64	131	60	129								
	SP	13	75	135	70	133								
	10 Petite	10												
	312 Petite	312												
	312 UZ Petite	312	M: 50 P: 60	M: 120 P: 123	M: 40 P: 50	M: 112 P: 117								
	312	312	SP: 70	SP: 128	SP: 60	SP: 122								
	312 UZ	312												
	13 UZ	13												

▲ Requires iCom ● Requires remote control option ■ Requires telecoil option

		Features																Phonak CROS											
		Binaural Directionality				Adaptive Intelligence				A Phonak has it all								CROS/ BiCROS											
		StereoZoom	autoZoomControl	ZoomControl	DuoPhone	UltraZoom Premium	UltraZoom Advanced	UltraZoom Standard	QuickSync	FlexControl	FlexVolume	SoundFlow Premium 4 programs	SoundFlow Advanced 3 programs	SoundFlow Standard 2 programs	User Preference Learning	SoundRecover	Real Ear Sound	WhistleBlock	NoiseBlock Premium	NoiseBlock Advanced	NoiseBlock Standard	WindBlock	EchoBlock	SoundRelax	Acoustically Optimized Vent (AOV)	Channels	Manual Programs	Compatibility	
Petite	Phonak Ambra																									20	5+2		
	Phonak Solana																										16	4+2	
	Phonak Cassia																										6	3+2	
microM microP SP	Phonak Ambra																										20	5+2	
	Phonak Solana																										16	4+2	
	Phonak Cassia																										6	3+2	
10 Petite 312 Petite 312 UZ Petite	Phonak Ambra																										20	5+2	
	Phonak Solana																										16	4+2	
	Phonak Cassia																										6	3+2	
312 312 UZ 13 UZ	Phonak Ambra																										20	5+2	
	Phonak Solana																										16	4+2	
	Phonak Cassia																										6	3+2	