

Phonak Baseo™ Q

Phonak Tao™ Q

The basics: great sound, great value





Great foundation

Over six decades, Phonak has established itself as the worldwide leading hearing health care provider with a solid reputation for innovation, advanced technology and reliability.

Through our global presence and network of local partners in over 100 countries we continue to reach out to enhance the quality of life for people with hearing loss.

With Phonak and your hearing care professional you are in the best hands on your journey to better hearing.

Roots in Switzerland

At Phonak we are proud of our Swiss heritage. At our headquarters near the Lake of Zurich, audiologists and engineers have designed, developed and produced industry-leading hearing aids since 1947. These include some of the smallest and most powerful hearing solutions available on the market. Phonak stands for perfection, highest quality, reliability and trustworthiness.



Today's modern hearing aids will surprise you. Gone are those bulky whistling products you may remember from the old days. Our newest generation of hearing aids are discreet, high-tech digital devices with excellent sound quality.

Even in our entry-level the hearing aids Phonak Baseo Q and Phonak Tao Q offer you the benefits of a state-of-the-art quality product:

- Great sound
- Comfortable listening in noise
- Proven reliability



Phonak Baseo Q

Phonak Baseo Q is available in four Behind-The-Ear models each with an attractive durable housing. Baseo Q covers a broad range of hearing needs from mild to severe hearing losses.

Discreet hair and skin colors



Phonak Tao Q

Phonak Tao Q comes in three models that are worn discreetly inside the ear. Each Tao Q is individually crafted to fit your own individual ear to ensure the best possible fit and comfort.



Please ask your hearing care professional about Tao Q colors and options.

Your benefits from Phonak Quest technology

Enjoy great sound

Phonak Baseo Q and Tao Q hearing aids are built on state-of-the-art Phonak Quest technology* also used in the most advanced Phonak products.

Baseo Q and Tao Q are computer programmed based on your individual hearing loss. Your hearing care professional can then further adjust the settings based on your individual preferences in 2, 3 or 4 frequency ranges. These adjustments are much more precise than any mechanically tuned hearing aid.

Enjoy comfortable listening – even in noise

Phonak Baseo Q and Tao Q in performance levels Q15 and Q10 include the single most effective technology in hearing aids to improve understanding in noisy situations: directional microphones. They reduce background noise and significantly improve your understanding of the conversation partner in front of you, for example in a busy restaurant.

Proven reliability

Baseo Q and Tao Q are based on our most robust designs – sold over a million times worldwide. And we improved them further to be your daily, highly reliable companion even when challenged by rain, sweat and dust.

Features comparison by performance level

Phonak Baseo Q15/Q10/Q5

Phonak Tao Q15/Q10

Key features	Q15	Q10	Q5
Frequency channels	4	3	2
Precision computer programming	•	•	•
Directional microphone option	•	•	-
Telecoil	•	-	-
Phonak BassBoost™	•	-	-

Other features

Phonak WhistleBlock	•	•	•
Phonak NoiseBlock Basic	•	•	•
Manual programs	4	3	2
Baseo Q color options	4	3	2
Tao Q color options	4	1	-
Roger compatibility	•	-	-
WaterResistant IP57**	•	•	•

**IP57 indicates that the hearing aid is water resistant and dust protected. It survived immersion in 1m of water for 30 minutes and 8 hours in a dust chamber as per the IEC60529 standard. IP57 apply to Phonak Baseo Q only.



* According to an international customer survey conducted by the Phonak Research Center, nine out of ten hearing care professionals would recommend Phonak Quest hearing aids.

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak.com

