


Introducing the new Phonak Virto™ B-Titanium



Even
more
discreet




- New Titanium FitGuide
- New module colors
- Smaller vent options
- Extra Retention
- Biometric Calibration

The smallest Phonak custom hearing aid just got even more discreet.

Ear impressions are static but ear canals are flexible. The new Titanium FitGuide takes this flexibility into account and indicates how deep the hearing aid can be placed in the ear canal.

Thanks to the new Titanium FitGuide, over **50%** of people get a deeper fitting Virto B-Titanium by an average of **2.5 mm**. That is the difference between discreet and super discreet.



What else is new?

New module colors



reddot award 2018
winner

Smaller vent options

With Virto B-Titanium we offer a new cleaning tool for smaller vents. It allows 0.8 mm vents to be built for optimal feedback thresholds, even for SP receivers.

Extra Retention

New optional shell surface technology.

Biometric Calibration

Precisely calibrated to the individual ear anatomy.