

Phonak Fast Facts



Phonak DECT cordless phone

What is it?

The Phonak DECT cordless phones are the ideal phones for use both at home or in a small office. Automatically, the phone signal is streamed wirelessly to both hearing aids. This provides unmatched sound quality and great speech understanding.

Apart from the streaming function, the phone also has a regular acoustic output and can therefore be used by non-hearing aid wearers too. Up to six handsets can be connected to one base which is particularly useful for small-office use. The phone features: One-click volume boost up to +15 dB, dedicated volume buttons, hands-free mode, a vibrate function and more. There are 2 version of the DECT cordless phone. DECT CP1 is used with clients who have Spice, Spice+ or Quest hearing aids. Phonak DECT II is for clients who have Phonak Venture and Phonak Belong™ hearing aids.



Why use it?

For hearing aid users, especially when they have a significant hearing loss, speaking on the phone can be challenging. Receiving the phone signal in both ears significantly improves speech intelligibility on the telephone, especially when background noise is present. Independent studies have shown that using the Phonak DECT phone, rather than a standard home phone, speech understanding can be improved by over 40%.¹

How is it used?

The DECT cordless phone is easy to use – no programming and no additional devices are required – the user can simply pick it up and talk! It can be used just like a normal phone – no need to manually change hearing aid programs. The user just holds it in the normal phone position (within 25 cm of the hearing aids) and the signal is streamed automatically.

¹ Appleton-Huber J. (2013) "Phonak DECT CP1 cordless phone", Field Study News. Phonak AG: 2013

For more detailed information visit www.phonakpro.com/dect