

When you can focus on work rather than on hearing, life is on

### Communicate

The modern workplace can be a challenging listening environment. It is also a place where successful communication is key in order to succeed.

### Participate

When affected by hearing loss, you become a master at finding ways around difficult hearing situations. This can be time consuming, emotionally demanding and often limits your ability to fully participate and engage with your workplace colleagues.

### Contribute

Roger wireless microphones can be used together with modern hearing aids, a solution that boosts your performance and makes it easier for you to interact at work.

### New groundbreaking MultiBeam Technology

Roger Select features MultiBeam Technology. By utilizing multiple microphones in six directions, speech from a 360° radius is calculated and compared. The direction with the best clarity is automatically selected.

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

[www.phonak.com](http://www.phonak.com)



## Roger™ wireless solutions at work

# Bridging the understanding gap

Even the most powerful hearing aids have limitations. The Roger portfolio consists of wireless microphones and receivers. They pick up the speaker's voice and transmit it to receivers connected to your hearing aids, so you can hear speech clearly over distance and in noise.

## Roger microphones



### Roger Table Mic II

Designed for meetings; several microphones can be connected to cover large meetings.



### Roger Select™

A versatile microphone which automatically selects the person who is talking and with Bluetooth® for phone calls.



### Roger Pen™

A handy microphone with portable design and Bluetooth for phone calls.

## Roger receivers



### Roger design-integrated

Click onto Phonak hearing aids and selected implants.



### Roger X

Universal miniature receiver.




### Roger MyLink

Neck loop receiver for hearing aids or implants with a telecoil.

Situations with background noise or where the sound source is at a distance are ones where Roger truly excels:

### Small meetings




Small meetings can be overwhelming for people wearing hearing aids. Place one microphone on the meeting table and hear your colleagues clearly.

### Large meetings



In large meetings multiple Roger Table Mic II will do an excellent job picking up all voices.


### Meetings with presenters



When the presenter wears a Roger microphone that works with the table microphone you can hear the speech clearly from all participants.


Additional situations where Roger can help overcome communication challenges:

### One-on-one conversation in noise




When surrounded by distracting noise hang a Roger microphone around your colleague's neck or clip it on to the collar.

### Multimedia



Video conferencing systems and Internet calls are frequently used at work. Simply connect your Roger microphone, using the audio cable.

### Phone calls



By connecting Roger Pen or Roger Select to your phone via Bluetooth you can hear phone calls directly through your hearing aids or implants.